



Beiwe Service Center Studies.beiwe.org Administrator Manual

Last updated: April 10, 2023

Beiwe Research Platform	3
Introduction	3
This Manual	3
Logging in/Changing Your Password	4
Setting Up Your Study	4
Surveys	4
Written Surveys	4
Absolute and Relative Survey Scheduling	8
Survey Branching	10
Special Characters in Surveys (Unicode)	15
Audio Surveys	16
Passive Data Collection	17
Allow Data Upload Over Cellular Plan	18
Consent via the Apps	18
Android	18
iOS	19
Checklist: Items to Provide to Onnela Lab	20
Adding Participants to Your Study	20
Downloading Data	21
Creating Your Data-Download Credentials	21
Download Data Tab	21
Download API	22
Frequently Asked Questions	22
Questions, Comments or Feedback?	22

Beiwe Research Platform

Introduction

Beiwe is a research platform designed to collect research-quality smartphone raw sensor and usage data. Our ultimate goal is to develop statistical and computational tools to enable the extraction of biomedical and clinical insights from smartphone data.

The Beiwe platform is comprised of a study portal, smartphone applications for Android and iOS operating systems, data storage, and data modeling and analysis tools designed for transparent, customizable, and reproducible biomedical research. The Beiwe research platform is used in both clinical and non-clinical research areas including psychiatry, neurology, surgery, oncology and basic neuroscience.

The Beiwe Research Platform is a project of the Onnela Lab and was built primarily by an NIH Director's New Innovator Award made to Dr. Jukka-Pekka "JP" Onnela.


This Manual

Welcome to the Beiwe Research Platform Administration Guide for Collaborators. This manual reviews how collaborators will utilize the Beiwe Research Platform's backend or "Study Portal" (<https://studies.beiwe.org>) to setup and review your study, add participants to your study and download your data. You will use studies.beiwe.org throughout the course of your research.

As a Beiwe user, you have the ability to do the following tasks through the <https://studies.beiwe.org> study portal website:


1. Create new Beiwe Patient IDs and temporary passwords to use to register your study participants
2. Create, edit, and deploy written surveys
3. Create, edit, and deploy audio surveys
4. Download study data

As illustrated in this diagram, this is the first step to get your study underway.



Logging in/Changing Your Password

Research collaborators will receive their login credentials from the Onnela Lab Research Administrator. This is a screenshot of the Login screen.



https://studies.beiwe.org

Log In

Username Password

(optional MFA code) Log in

Please immediately change your password to a secure password of your choice by clicking on the **“Manage Credentials”** tab. Complete the **“Reset Your Password”** section and click on the **“Change Password”** button.

Reset Your Password

Don't forget your password!

Beiwe does not store your email address, and has no other means of getting in contact with you.

If you forget your password you will have to contact a Beiwe administrator to have it reset.

The primary Beiwe administrator is Kenzie Carlson, she can be reached at kcarlson@hsph.harvard.edu.

Passwords require all of the following:

- at least one lower case letter
- at least one upper case letter
- at least one numeric character
- at least one special character (! @ # \$ % ^ [] etc.)

Current Password:

New Password:

Confirm New Password:

[Change Password](#)

Setting Up Your Study

Surveys

Written Surveys

Choose the name of your study from the list on the “**Manage Studies**” screen or click “**View Study**” in the top left side of the screen. Scroll down past the Participants section to the “**Surveys**” section and click “**Create New Survey**” to create a survey. If you have multiple surveys for a study, they will appear in a list.

Surveys

[Create New Survey](#)

This study does not have any Surveys.

Audio Surveys

[Create New Audio Survey](#)

This study does not have any Audio Surveys.

Edit Survey

Test_1

The survey name is not visible to the study participant.

Schedule

☐ Show survey notification immediately when the user downloads this survey. (Will also trigger this survey right after registration.)

Schedule type:

☒ Weekly
 ☐ Relative
 ☐ Absolute

Weekly schedules send surveys at the same time every week

Add new scheduled time (type numbers into the pop-up time picker):

Currently scheduled times:

Sunday

Monday

Tuesday

Wednesday

Thursday

Friday

Saturday


Questions

Your changes will not be saved until you click "Save and Deploy."

☐ **Randomize** which questions are displayed each time the participant takes a survey. This also randomizes the order in which the questions are displayed. It also causes the survey to ignore all question display logic so that all questions have an equal chance of appearing.

☐ **Make tracking survey always available** if checked, the survey will be available to the participants at all times regardless of schedule

To add a question to the survey, click the “**Add Question**” button at the bottom of the page.



To create and edit questions, fill in the appropriate fields in the pop-up, including Question Text, Question Type, and Display Logic. Then, click the “**Add Question**” button.

- i. Question Types include:
 1. Informational Text
 2. Numerical Slider
 3. Radio Buttons
 4. Checkboxes
 5. Free Text Response

After making any changes, click the “Save and Deploy” button to save the changes.

Edit the timing of the surveys by using the scheduling system under the **Schedule** section. Click the “Add Time” button to set the time.

Edit Survey

Test_1

The survey name is not visible to the study participant.

Schedule

☐ **Show survey notification immediately** when the user downloads this survey. (Will also trigger this survey right after registration.)

Schedule type:

Weekly **0** Relative **0** Absolute **0**

Weekly schedules send surveys at the same time every week.

Add new scheduled time (type numbers into the pop-up time picker):

Currently scheduled times:

Day	9:00 AM	5:00 PM
Sunday	<input type="button" value="X"/>	<input type="button" value="X"/>
Monday	<input type="button" value="X"/>	<input type="button" value="X"/>
Tuesday	<input type="button" value="X"/>	<input type="button" value="X"/>
Wednesday	<input type="button" value="X"/>	<input type="button" value="X"/>
Thursday	<input type="button" value="X"/>	<input type="button" value="X"/>
Friday	<input type="button" value="X"/>	<input type="button" value="X"/>
Saturday	<input type="button" value="X"/>	<input type="button" value="X"/>

After any change, click the “Save and Deploy” button to save the changes.

Schedule

☐ **Show survey notification immediately** when the user downloads this survey. (Will also trigger this survey right after registration.)

Schedule type:

Weekly **14** Relative **0** Absolute **0**

Weekly schedules send surveys at the same time every week.

Add new scheduled time (type numbers into the pop-up time picker):

Currently scheduled times:

Day	9:00 AM	5:00 PM
Sunday	<input type="button" value="X"/>	<input type="button" value="X"/>
Monday	<input type="button" value="X"/>	<input type="button" value="X"/>
Tuesday	<input type="button" value="X"/>	<input type="button" value="X"/>
Wednesday	<input type="button" value="X"/>	<input type="button" value="X"/>
Thursday	<input type="button" value="X"/>	<input type="button" value="X"/>
Friday	<input type="button" value="X"/>	<input type="button" value="X"/>
Saturday	<input type="button" value="X"/>	<input type="button" value="X"/>

Questions

Your changes will not be saved until you click "Save and Deploy"

- If you would like the survey to appear when the participant in the study first downloads the app, click the box “**Show survey notification immediately** when the user downloads this survey. (Will also trigger this survey right after registration.)” If no days/times are scheduled

for the survey it will only show once – when the participant first enrolls. (However, if the participant gets a new phone or needs to delete/reinstall Beiwe, if they use the same patient ID, they will see this survey again. You may/may not want to instruct participants to complete it or just skip through the questions and submit so as to delete it.)

Schedule

☒ Show survey notification immediately when the user downloads this survey. (Will also trigger this survey right after registration.)

Schedule type:

Weekly 14

Relative 0

Absolute 0

- If you would like to use the randomized survey functionality, click on the “**Randomize**” or the “**Randomize without replacement**” button on the bottom of the page and specify how many questions you’d like to be included in each survey. “**Randomize without replacement**” means all questions will be presented before starting over, while “**Randomize**” will just choose the specified number of questions each time the survey is presented.

+ Add Question

☒ **Randomize** which questions are displayed each time the participant takes a survey. This also randomizes the order in which the questions are displayed. It also causes the survey to ignore all question display logic so that all questions have an equal chance of appearing.

☐ **Randomize without replacement**, i.e., on each survey, only display questions that haven’t appeared yet on a previous survey, until all questions have appeared in a survey, and then start over

How many questions are displayed in each survey, out of the 5 questions currently in the question bank for this survey

☐ **Make tracking survey always available** if checked, the survey will be available to the participants at all times regardless of schedule

Save & Deploy

Delete survey

Absolute and Relative Survey Scheduling

New survey scheduling features are available on the new versions of the Beiwe app for both Android and iOS. These features include relative survey scheduling, absolute survey scheduling, and setting the time zone of a study.

Relative Survey

Relative schedules send surveys a specified number of days before or after a date and time that is specific to each participant, which is known as an **intervention**. Surveys arrive at the designated date and time in the **designated study-wide time zone**.

1. Create an intervention
 - a. On the study’s page, click the “**Edit Interventions**” button under **Study Configuration** to navigate to the Interventions page.

Surveys

[Create New Survey](#)

Test_2

Survey ID #kzvS1NE3zEUz3Wxy8VjK71eY
Updated: 2022-03-24 7:34AM (EDT)

Test_1

Survey ID #Pf8RVZ1TqWmx1V9yVBzO8F7G
Updated: 2022-03-24 7:25AM (EDT)

Audio Surveys

[Create New Audio Survey](#)

This study does not have any Audio Surveys.

Image Surveys

[Create New Image Survey](#)

This study does not have any Image Surveys.

Study Configuration

[Edit Interventions](#)

Configure Interventions for use with Relative survey schedules

[Edit App Settings](#)

Configure types and quantity of the passive data streams this study collects, and the wording that study participants see in the app

[Edit Custom Fields](#)

Edit custom tags for organizing participants in your study.

b. Create a new intervention by typing in a name and clicking “**Add New Intervention**”.

Edit which Interventions are on this study

Interventions	Edit Interventions	Remove Interventions
Intervention_1	Edit Intervention	Remove Intervention

Add new Intervention

Intervention Name:

[+ Add New Intervention](#)

c. Back on the study page, the new intervention will appear as a column in the **Participants** section.

Participants

Total participants ever registered on this study: 0

Show entries

Search:

Creation Date	Patient ID	Phone registered	Phone OS	Intervention_1
2022-03-24	vbub3sbg	false		
2022-03-24	vbgj1rpr	false		
2022-03-24	aolh2rs1	false		

d. Click on each Patient ID to set the value of the intervention for each participant.

Activity Information

Phone operating system: Unknown?

Beiwe App Version:

Time of Registration:

Time of last check-in for updated survey content:

Time of last file upload event:

Time of last notification credential update:

Time of last device settings update:

Survey Notifications

Total notifications or attempts: 0 [Notification History](#)

Enable Easy Enrollment

Easy enrollment may be enabled on a per-participant basis. When easy enrollment is enabled for a study or a participant the Temporary Password that is otherwise required at registration is ignored. Any value that passes any validation performed by the app is allowed.

Generate a new password for this participant.

At registration the participant enters their current password into the temporary password field, and then sets a new password that is used to unlock the app and access new surveys, etc. The participant may change this password from their device at any time. You cannot view the participant's current password.

Un-Register Device

There is no device registered for this participant

Resetting a participant's registration password is a two-step process. You should only reset a participant password if you know there is no device registered and uploading data. Click this button if you need to register a new participant or reregister an old participant who has forgotten their password.

Permanently Retire Participant

This action retires this participant entirely, blocking any further data from being uploaded, regardless of the source device. This does not affect data already uploaded. Attempts to register a new device with this participant will fail. If data is uploaded for this participant the device will be told to delete that data.
This operation is final and cannot be undone.

Interventions

Intervention: mm/dd/yyyy @

2. Create a new survey. On the survey page, navigate to the **“Relative”** tab. Create a timing for the absolute survey using the options on that tab, relative to the selected intervention. Options include “Day of”, “Days After” and “Days Before”.

Schedule

☐ **Show survey notification immediately** when the user downloads this survey. (Will also trigger this survey right after registration.)

Schedule type:

Weekly 14 Relative 0 Absolute 0

Relative schedules send surveys a specified number of days before or after an intervention date that is specific to each participant

Add new scheduled time followed by how many days before or after the intervention event (type numbers into the pop-up time picker):

08:00 AM 🕒 ✓ Days After Days Before Day Of intervention_1 + Add time

Note: Scheduled events based on Relative schedules occur is in the past (usually within 5 minutes). Beiwe will not send survey notifications if there was already one sent for the date and time of the trigger.

Currently scheduled times:

Absolute Survey

Absolute schedules send out surveys to participants at an exact date and time. Surveys arrive at the designated date and time in the **designated study-wide time zone**. If a patient registers after the date of the absolute survey, it will be delivered upon registration.

1. Create a new survey. On the survey page, navigate to the **“Absolute”** tab. Create a timing for the absolute survey using the options on that tab.

Schedule

☐ Show survey notification immediately when the user downloads this survey. (Will also trigger this survey right after registration.)

Schedule type:

Weekly 14 Relative 0 Absolute 0

Absolute schedules send out surveys to participants on a specified date

Add new scheduled time followed by the scheduled date (type numbers into the pop-up time picker):

08:00 AM 03/24/2022 +Add time


Note: Absolute schedules created with a time in the past will trigger immediately (usually within 5 minutes). Beiwe will not send survey notifications if there was already one sent for the date and time of the schedule.

Currently scheduled times


Set the time zone of a survey

Set the time zone for a survey, on which all surveys will be sent out. **Note:** It is important to set the time zone at the beginning of the study and do not change it after sending surveys to participants, as this will cause surveys to be reset.

1. From the study page, click **“Edit Study”** in the upper left side of the screen.



2. Scroll down until you see the option to **“Set Study Timezone.”** Select the appropriate time zone and click **“Update Timezone”**.




Survey Branching

Any questions that you would like to use to determine branching logic must have a definitive answer, thus must use the following question types:

- slider questions,
- radio button questions, and
- numeric open response questions.

Beiwe's remaining question types, such as information or checkbox, do not result in a unique outcome thus cannot be used for branching decisions. However, they can be used in the normal flow of your surveys without being referred to later. Set up your questions as you always have:

- Enter your question text
- Choose a question type (as noted above, depending on how you're using the question)
- Add your answer choices.
- Press **"Add Question"** to save



Q1

Always Displayed

In the last 7 days did you experience any abdominal pain?


Question Type: Radio Button

Options: Never , Rarely , Occasionally , Frequently , Almost constantly

Edit Delete Move Up Move Down

Questions in surveys are now designated as "Always Displayed" or "Conditionally Displayed."

As you add your answer choices to a radio button survey question, note that a number is now noted to the left of your list of the choices, as shown in the screenshot below. If you'll be using this question for survey branching, take note of these numbers to use later on when you set up your branching logic.



Edit Question

Question Text: In the last 7 days, how much did your Abdominal Pain interfere with your usual or daily activities?

Question Type: Radio Button

Answer Choices: Add new answer choice

0: Not at all x

1: A little bit x

2: Somewhat x

3: Quite a bit x

4: Very much x

If you're using a slider or numeric open response question type, the specific number entered on either can be used for your subsequent conditional questions.

“Always Displayed” and “Conditionally Displayed” questions can be used as the basis for survey branching. To set up the conditions for displaying a “Conditionally Displayed” question, the following four operator blocks are available:

- Add “NOT” block
- Add “OR” block
- Add “AND” block
- Add Conditional

One or more of the operators can be configured to implement your branching logic. Click on one or more of these new buttons in the **Edit Question** screen to get started. “**Add Conditional**” is a single operator as noted when used and as can be seen in the following screenshot.

Edit Question

Question Text:

Question Type: Radio Button ▾

Answer Choices: Add new answer choice

0:

1:

2:

3:

4:

Display this question if:

You will be limited to a single conditional statement if you do not have an AND or OR block as the first element.

Answer to Q2 ▾ > ▾ ×

The NOT, OR and AND blocks can be configured with multiple conditional statements and can be nested within each other. See an example in the following screenshot. Remember to click the “**Add Question**” button to save your questions and the “**Save and Deploy**” button to save the survey when you’re done and as you’re adding new questions. Add your survey delivery times as normal before your final **Save and Deploy**.

[illegible]

Studies.beiwe.org Administrator Manual

Questions

Q1
Always Displayed

Did you go on a run today?
Question Type: Radio Button
Options: No , Yes

Edit Delete Move Up Move Down

Q2
Conditionally Displayed

How many miles did you run?
Question Type: Slider
Minimum: 0
Maximum: 13

Edit Delete Move Up Move Down

Q3
Conditionally Displayed

Get off the couch!
Question Type: Informational Text

Edit Delete Move Up Move Down

Edit Question

Question Text:

How many miles did you run?

Question Type:

Slider

Minimum Value:

0

Maximum Value:

13

Display this question if:

You will be limited to a single conditional statement if you do not have an AND or OR block as the first element.

Answer to

Q1

==

1

When will I get an error when entering branching logic?

You won't get errors when you save contradictory logic, like the expression " $((Q1 == 0) \text{ AND } (Q1 == 1))$ ". You will have to correct logic errors like this on your own when you notice a survey is not giving you the results you're expecting. Errors only appear when Beiwe can't parse the logic, for example, if you're referencing a question that doesn't exist anymore (e.g., if you referenced Q5 and then deleted it, even if there's a new Q5) or if you're referencing a question that doesn't have a numeric answer (like a checkbox question or alphanumeric response). If the logical expression can be parsed, even if the logical expression always evaluates to *False*, it won't show any errors when you click **"Save and Deploy"**.

What happens to questions that aren't displayed?

If any comparative operator (<, <=, ==, >=, or >) references a question that wasn't answered, the expression should evaluate as False. For example, if you're evaluating (Q3 < 1) AND (4 > 3) and Q3 wasn't answered, it should evaluate as (False) AND (4 > 3). If you're evaluating NOT(Q3 == 1) and Q3 wasn't answered, that should evaluate as NOT(False), which equals True.

If Q3 never displayed, it cannot have been answered. If Q3 displayed, and the participant answered it, and then the participant went back in the survey and changed a previous answer so that Q3 shouldn't have displayed, the answer to Q3 is still there to be operated on. So if you want Q4 to never display if Q3 didn't display either, you could add a copy of the Q3 display logic to Q4's display logic inside an AND block.

Entering long text responses: How to write responses so the app shows a new line.


You need to use two newlines when entering your text in the “**Question Text**” box in order to get a paragraph break to display in the app. A single newline is ignored.

Beiwe Randomize feature

Beiwe's randomize feature does not work with survey branching. If you click the “**Randomize**” box when configuring your survey, it will disable any branching questions in the survey.

Special Characters in Surveys (Unicode)

Beiwe allows for unicode characters to be added to questions to allow for surveys in other languages. However, at this time the hardcoded language in Beiwe is only available in English. See an example of a deployed survey below on an iPhone XS (OS 14.0.1):



Audio Surveys

Audio surveys are created and scheduled in the same way that surveys are scheduled.

The researcher can enter the **“Voice Recording Prompt”** directly into the **Edit Survey** section. This will appear on the screen when the participant is directed to the audio survey section from the app notification. Wording for the **“Voice Recording Prompt”** is completely customizable. **Researchers should hit the “Save and Deploy” button after any edits or additions.**

Voice Recording Prompt

Technical Settings

Audio recordings contain a single channel and are made using the active microphone on the device, which may be on a connected headset. Older versions of the Beive App will make compressed recordings at 64Kbps.

Compressed audio files use AAC compression. 64Kbps should be sufficient for a vocal recording under nearly all situations, and produces roughly **one half of a megabyte** of data for each minute of audio. The sample rate of a compressed recording is always 44,100Hz.

Uncompressed audio recordings produce wav files. A wav file with a sample rate of 44,100 ("CD quality") has a bit rate equivalence of 706Kbps, or roughly **five megabytes** of data for each minute of audio.

☒ Compressed audio ☐ Uncompressed audio 64Kbps ▼

☐ Make audio survey always available if checked, the survey will be available to the participants at all times regardless of schedule

Save & Deploy Delete survey

Technical Settings: Researchers should use the **“Compressed audio”** feature as a default, unless their study has specific needs related to uncompressed audio.

Uncompressed audio recordings sample at a higher rate (16,000, 22,000 or 44,100 Hz) and can reach the quality equivalent to CD quality. However, the cost is excessive data storage. Roughly 5MB of data would be collected for each minute of uncompressed audio at the 44,100 Hz level.

Although uncompressed audio takes up quite a lot of space, compression leads to compromised audio features (even if not audible to a person), and this may have data implications depending on your study aims. Please consult with the Onnela Lab if you would like to discuss this further.

Technical Settings

Audio recordings contain a single channel and are made using the active microphone on the device, which may be on a connected headset. Older versions of the Beive App will make compressed recordings at 64Kbps.

Compressed audio files use AAC compression. 64Kbps should be sufficient for a vocal recording under nearly all situations, and produces roughly **one half of a megabyte** of data for each minute of audio. The sample rate of a compressed recording is always 44,100Hz.

Uncompressed audio recordings produce wav files. A wav file with a sample rate of 44,100 ("CD quality") has a bit rate equivalence of 706Kbps, or roughly **five megabytes** of data for each minute of audio.

☒ Compressed audio ☐ Uncompressed audio 64Kbps ▼

☐ Make audio survey always available if checked, the survey will be available to the participants at all times regardless of schedule

Passive Data Collection

Some settings, including the Passive Data settings, and the consent text, can only be edited by an Onnela Lab member but you will be able to view the settings. Since the Beiwe Research Platform is geared toward collection of research-grade passive sensor data, the app itself is meant to be minimalistic.

Below is a screenshot of the configurable passive data settings supported by the platform. These are the current default settings for a study. However, this may change at any time and you may request specific settings to support your study aims. When considering the passive data settings and tradeoffs, it may be helpful to note that the global positioning system (GPS) sensor setting is more of a battery issue in that this sensor utilizes the phone's battery capacity while the Accelerometer produces a lot of data, so would impact data storage capacity.

HSPH_Kennedy_Test_Study / Study Settings / App Settings

Data Stream Settings

Accelerometer, GPS, Gyro, and Magnetometer are available on Android and iOS.
Power State is present on both platforms, but the available data differs.

Device Motion, Reachability, and Proximity are iOS-only.
Ambient Audio, Wifi, Calls, Texts, and Bluetooth are Android-only.

GPS:	<input checked="" type="checkbox"/>	GPS on duration:	<input type="text" value="60"/>	GPS off duration:	<input type="text" value="600"/>
Magnetometer:	<input type="checkbox"/>	Magnetometer on duration:	<input type="text" value="60"/>	Magnetometer off duration:	<input type="text" value="600"/>
Device Motion:	<input type="checkbox"/>	Device Motion on duration:	<input type="text" value="60"/>	Device Motion off duration:	<input type="text" value="600"/>
WiFi:	<input checked="" type="checkbox"/>	WiFi Log Frequency:	<input type="text" value="300"/>		
Accelerometer:	<input checked="" type="checkbox"/>	Accelerometer on duration:	<input type="text" value="10"/>	Accelerometer off duration:	<input type="text" value="10"/>
				Accelerometer Frequency:	<input type="text" value="10"/>
Gyro:	<input type="checkbox"/>	Gyro on duration:	<input type="text" value="60"/>	Gyro off duration:	<input type="text" value="600"/>
				Gyro Frequency:	<input type="text" value="10"/>
Ambient Audio:	<input type="checkbox"/>	Ambient Audio on duration:	<input type="text" value="600"/>	Ambient Audio off duration:	<input type="text" value="600"/>
		Ambient Audio Bitrate:	<input type="text" value="24Kbps"/>	Ambient Audio sampling rate:	<input type="text" value="44,100"/>
Power state:	<input checked="" type="checkbox"/>	Proximity:	<input type="checkbox"/>	Reachability:	<input checked="" type="checkbox"/>
				Calls:	<input checked="" type="checkbox"/>
				Texts:	<input checked="" type="checkbox"/>
<i>The Bluetooth data stream is synchronized across all devices in your study. It scans at a specific time within a defined period. An On duration of 5 minutes with a Total duration of 1 hour and a global offset of 15 minutes will trigger a bluetooth scan for 5 minutes once an hour, every hour, at 15 minutes past the hour.</i>					
Bluetooth:	<input type="checkbox"/>	BT on duration:	<input type="text" value="60"/>	BT total duration:	<input type="text" value="300"/>
				BT global offset:	<input type="text" value="0"/>

A participant in a study with all data streams turned on could generate about 1 gigabyte of data per month. Most of this is accelerometer data, since every accelerometer datapoint is about 50 bytes of data, and the accelerometer records about 5 to 100 datapoints per second. Other data streams produce less data simply because they record less frequently. This would likely be a worst-case scenario as more recently we've been limiting the accelerometer data to reduce these numbers.

The Android and Apple iOS phones support different passive data collection items. The table below shows which data streams can be collected on each phone. The Beiwe study settings will be ignored if the phone does not support an enabled setting. Thus, if **"Phone Calls"** is enabled, this data will only be collected from Android phones since Apple does not permit accessing this data.

Supported Data Streams

Sensor/Data Stream	Android	iOS
Accelerometer	Yes	Yes
Android Log File	Yes	No
Bluetooth	Yes	No
Device Motion	No	Yes
GPS	Yes	Yes
Gyroscope	Yes	Yes
Identifiers	Yes	Yes
iOS Log File	No	Yes
Magnetometer	No	Yes
Phone Calls	Yes	No
Power State	Yes	Yes
Proximity	No	Yes
Reachability	No	Yes
Survey (Audio)	Yes	Yes
Survey (Text)	Yes	Yes
Text Messages	Yes	No
Wi-Fi	Yes	No

Allow Data Upload Over Cellular Plan


Beiwe typically relies on a reliable Wi-Fi connection to upload data collected from a study participant. However, for studies that would like the data more immediately, or have other specific study aims, a study can be configured for data upload to go over the participant's cellular plan. As this can be costly, collaborators generally reimburse study participants for the cost of their data plans or only accept participants with unlimited data plans. Study participants will have different providers for their cell phone plans and each provider offers different terms for their data plans and tends to change the terms frequently.

Consent via the Apps

While most collaborators choose to consent study participants using the traditional paper-based methods, Beiwe can also be used to provide consent electronically if approved by your IRB. The Android and iPhone have two distinct consent processes. You will need to work with Onnela Lab to set up both to accommodate your study participants. However, the text can be re-used.

Android

The default text for the Android Consent Form is shown on the screenshot below. The Researcher will go through the Required Readings (detailed below) with the participant, after which the participant will hit the “**I Consent**” button, and the app will be activated on the phone.



iOS

The iOS Beiwe app offers seven headings in the consent process of which you may use all, some or none. The default set up is that all fields are blank and that none of these screens are shown during registration. However, if you would like to use these screens for sharing information or gaining consent, just provide your desired text for one or more of these topics to Onnela Lab.

The seven optional topics are:

- Welcome
- Data Gathering
- Privacy
- Time Commitment
- Study Survey
- Study Tasks
- Withdrawing

Each topic has two screens. The first screen is designed to give a one or two sentence explanation as shown above and the second screen provides the “learn more” information.

After stepping through each of the topics that are populated with text, the “learn more” text is compiled into one document as shown below and is presented to the participant for review and to ask for their consent to join the study.

Sample Consent Screen Review & Acceptance

Checklist: Items to Provide to Onnela Lab

After the IRB is finalized, the study can get underway. Collaborators should meet with the Onnela Lab to ensure the Passive Data Settings are agreed upon and the following items should be provided to the Onnela Lab:

- ☐ Title for Project
- ☐ PI should send team's names and email addresses that will need access to the Beiwu study
- ☐ Copy of the IRB approval
- ☐ Copy of the IRB materials
- ☐ Text for consent – Android, if applicable
- ☐ Text for consent – Apple, if applicable
- ☐ Text for the app's "Call My Clinician" button
- ☐ Approximate number of patients/study participants anticipated in study
- ☐ Approximate duration of the study
- ☐ "Allow Data Upload Over Cellular Plan" Setting

Adding Participants to Your Study

Researchers will enroll new participants in the study by creating and assigning Beiwu Patient IDs.

Click the **"Add New Participant"** button at the top of the screen, and a new 8-character Beiwu ID and temporary password will be generated on the screen. You can also Click the **"Add Many New Participants"** button on the study screen, and a window will pop up asking the number of patient IDs you would like to generate and the name of the file you would like the IDs to export to.

The new Beiwu ID and temporary password should be given to a study participant to be used when they register for your study from the Beiwu app on their phone.

iPhone Beiwe registration screen

Android Beiwe registration screen

It is the researcher’s responsibility to keep track of the Beiwe Patient ID and the participant’s identity. In keeping with IRB requirements, the Beiwe app ties all data collected to the Beiwe Participant ID and does not know the participant name or other identifying information.

Downloading Data

Creating Your Data-Download Credentials


The first time you want to download data from studies.beiwe.org, researchers will need to go to the **“Manage Credentials”** page, and under the **“Generate New Data-Download API Access Credentials”** heading, **click the yellow button (“Reset Data-Download API Access Credentials”).**

Keep these credentials in a secure/password protected area. You will need these credentials to access data each time you want to download data; these will be your **same credentials to access data for all the studies you have access to in the studies.beiwe.org system.** You only need to generate the data download credentials the first time you access data (you do NOT need to generate new credentials each time you access data). Once you have created and saved your credentials in the above step, for all future data download needs you can go directly to the step below and use your saved credentials.

Download Data Tab

Go to the “**Download Data**” page and enter your access key and secret key and select the study, patient ID, data types, and date ranges desired. Use CtrlA (control all) to select all data types for a particular patient ID. The End Date and Time will default to the current date if not entered.

Here’s an example of what the screen will look like.

The screenshot shows a web form for downloading data. It includes fields for 'Access Key' and 'Secret Key'. Below these are three columns: 'Study' with a dropdown menu showing 'HSPH_Test study_Maria'; 'Patients' with a list of IDs including '2khw7ost', '47il', 'brxyvprf', 'm8izaegl', 'mhqgq99', 'ooazoro5', 'pb85sto2', and 'ws4zyaov'; and 'Data Types' with a list including 'accelerometer', 'bluetooth', 'calls', 'gps', 'identifiers', 'app_log', 'power_state', 'survey_answers', 'survey_timings', 'texts', 'audio_recordings', and 'wifi'. At the bottom, there are date and time pickers for 'Start Date and Time' (10/19/2016 12:00 AM) and 'End Date and Time' (10/21/2016 12:00 AM), followed by a green 'Download Data' button.

Data Download Screen Example

Click “**Download Data**” and a zip of the requested csv files will be created. After downloading data, you will need to **refresh the page and re-enter your access and secret keys to download more data**. This is how the web page was designed for security purposes.

The format of the data files downloaded will vary by sensor according to the specific data being collected, however, the same metadata is provided for all raw data collected or any event recorded. Specifically, the first two columns for all raw data will provide the Unix timestamp and the UTC time to aid with data analysis. For more information, including an explanation of the sensors and data collected by Beiwe, please see the Beiwe Data Privacy and Security manual.

Download API

A Beiwe API is also available to download data directly to collaborators’ servers and databases. We recommend the API be used for downloading large amounts of data. Please contact Onnela Lab for more information.

Frequently Asked Questions

Please visit: <https://github.com/onnela-lab/beiwe/wiki/FAQs-and-Beiwe-Troubleshooting>

Questions, Comments or Feedback?

If you have any questions comments or feedback, please reach out to the Beiwe administrator: Kenzie Carlson, Research Manager: kcarlson@hsph.harvard.edu